



medibleblue™

**Surgical Type-IIR
Medical Masks**

PACKAGING AND SAFETY Product Specification Sheet

Product Codes:

MBSFM3

Sizes:

17.5 x 9.5 cm x 19.5

Features & Benefits:

- ✓ Synthetic Blood Penetration Test
- ✓ Microbial Cleanliness
- ✓ Differential Pressure
- ✓ Standards

- Tested in Controlled environment for better efficiency.
- Medical Type-IIR
- Suitable for jobs which requires maximum protection and grip
- 100% Ignited but extinguished results
- Comfortable, Easy fit and provide excellent protection and fit
- Excellent comfort, exact fitting and sensitivity with ear-loop
- Tested Against ISO and EN standards

Product Specifications:

PASSED

Medical TypeIIR

Physical Properties:

(typical):

Test Item	Test Result
Bacterial Filtration Efficiency (BFE) Test	≥98%
Fluid Pressure (mmHg)	Passed
Microbial Cleanliness	Passed
synthetic Blood Penetration	Passed

Quality Standards:

Conforms to standards and Manufactured in accordance with

- ✓ EN1468:2019 +AC:2019(E) Annex B
- ✓ EN1468:2019 +AC:2019(E) Annex C
- ✓ EN1468:2019 +AC:2019(E) Annex D
- ✓ EN 14683:2019+AC:2019(E)
- ✓ TUV reports
- ✓ ISO 22609:2004
- ✓ EN ISO 11737-1:2018
- ✓ SGS report
- ✓ CE
- ✓ EN ISO 14971-2019
- ✓ ISO 10993-1:2018
- ✓ EN ISO 15223-1:2016
- ✓ EN1041 : 2008



We proudly present you for your needs **Medical TypeIIR- Medical Mask**
Fulfilling your standard needs and worry-less work attitude to get the job done.
With the nose fit comfort and Ear-loop ease, it brings more to the table when it comes to Performance.

PRESENTED BY : **Myst Enterprise Ltd**
Christchurch-New Zealand

CONTACT US: +64 27 945 6289
info@mystenterprise.com
www.mystenterprise.com



- ✓ For Sensitive Skin
- ✓ No Stimulus to Human Skin
- ✓ Synthetic Blood Penetration Test Passed
- ✓ Nose Bar Adjust
- ✓ Easy Breathability, Perfect Fitting, Soft and comfortable

ATTENTION

1. Ensure the packaging is complete before using. Do not use if the packaging is damaged.
2. Use the respirator within the validity period.
3. Disposable respirator, do not reuse.
4. In case of allergy, stop using immediately.
5. Discard the respirator according to the regulations of environmental protection department.



**DISPOSABLE
MEDICAL
FACEMASK**



MBSFM3
50 Pcs

**Inert on Human Skin
Provide Excellent Protection**

Model C
Type IIR

For Sensitive Skin

BFE ≥ 98% 100% Kiwi Owned & Operated